



HOTLINE® Fluid Warmer


THE TOP THREE REASONS YOU USE THE LEVEL 1® HOTLINE® FLUID WARMERS

1

#1 in Hypothermia Prevention

HOTLINE® delivers normothermic blood & fluids all the way to the patient so your hospital avoids unnecessary complications and costs.

Blood and fluids delivered at normothermic temperatures without patient line cool down.




2

Outperforms Other Fluid Warmers

The HOTLINE® triple lumen tubing eliminates patient line cool down at routine flow rates, keeping your patients warmer and your PACU less crowded.

Outperforms hotplate technologies at routine flow rates.




3

You Get the Most Out of Every Infusion

The low priming volume means more fluid, blood and medications are delivered directly to the patient, minimizing waste.

The HOTLINE® patented triple lumen patient line maintains a layer of 42°C circulating solution around IV line.



Now we offer three new reasons to enhance your existing HOTLINE® Fluid Warmer, because the best fluid warmer just got better



Get HOTLINE® 3 for Free

BY CONVERTING YOUR EXISTING HOTLINE® FLUID WARMERS

Contact your Smiths Medical sales representative or customer service to schedule your "HOTLINE® 3 for Free" upgrade.

For more information visit our website:

www.smiths-medical.com/hotline3



ORDERING INFORMATION

DESCRIPTION	Part Number
HOTLINE® 3 Fluid Warmer	HL-390
HOTLINE® 3 Fluid Warming Set 20/case	L-370
HOTLINE® 3 Temperature Check	HLTA-390



1 Easier Set Up

The new disposable set (L-370) connects into the new mounting block with an audible click to let you know it's securely in place. It can now be disconnected quickly, easily and with one hand.

2 Easier Portability

The new disposable set helps keep recirculating solution in the tubing when disconnected from the warmer, so that you can transfer the tubing with your patient from the OR to the PACU and connect into another HOTLINE® 3 Fluid Warmer to continue patient warming.



3 Easier Maintenance

The HOTLINE® 3 conversion eliminates your monthly preventative maintenance requirement*, since O-rings no longer require lubrication.

*When using the Hydrogen Peroxide recirculating solution protocol.

Smiths Medical ASD, Inc.
Norwell, MA 02061, USA
Phone: 1-214-618-0218
Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Smiths Medical Canada Ltd.
Markham, Ontario, Canada, L3R 4Y8
Phone: 905-477-2000
Toll-Free: 1-800-387-4346

The Smiths Medical and Level 1 design marks and HOTLINE are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

Smiths Medical is part of the global technology business Smiths Group plc.
©2010 Smiths Medical family of companies. All rights reserved. TM194831EN 4/10

